



Single Chamber 8 Channels



Double Chamber 16 (8*2) Channels

POCtia

Poet Single Serving Automatic Chemiluminescence Analyzer

Accurate

- Enzymic (ALP) chemiluminescence, accurate full quantification
- Quality control and calibration consistent with large luminometers
- $CV \leq 5\%$
- No hydraulic circuit, ultra-low failure rate

Fast

- Result within 15 min
- 32(single chamber)/64(double chamber)tests per hour
- 8 channels/chamber in parallel, quickly process a large number of specimens
- Each chamber operates independently and flexibly

Intelligent

- Humanized operation interface, support two-way LIS connection
- Automatic item identification, saving manpower
- Tips、 magnetic sleeve grab and kick detection

Economic

- Single package without reagent waste
- Ready to disassemble and use, no reagent loss
- No need for additional consumables such as cuvettes、 Tips
- LCD touch screen, built-in computer, printer

Test Menu

Cardiac Marker

- hs-cTnI
- NT-proBNP
- D-DIMER
- CK-MB
- MYO



Inflammation

- IL6
- PCT
- CRP
- SAA




Nutritional

- FER
- FA
- VB12
- 25-OH-VD



Tumor Marker

- CEA
- AFP
- CA125
- CA19-9
- CA15-5
- T-PSA
- F-PSA



Fertility

- Total-hCGp
- PROG
- FSH
- PRL
- LH
- E2
- TESTO
- AMH



Thyroid

- **TSH**
- TT3
- TT4
- FT3
- FT4



Diabetes

- HbA1c



Allergy

- IgE



AFRICAN ALLIANCE MEDICAL

Tel: +237 675 322 450

Website: www.afralmed.com

Location: Rue Nsimeyong, Yaoundé

E-mail: contact.afralmed.org